

蛍光試薬	ロングパス	バンドパス
	シングルチャンネルモデル名	シングルチャンネルモデル名
Acridine Orange + DNA	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Acridine Orange + DNA and RNA	SXJ-5822FLS-19002	
Alexa Fluor 350™	SXJ-5822FLS-19000	SXJ-5822FLS-39000
Alexa Fluor 488™	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Alexa Fluor 546™	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Alexa Fluor 555™	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Alexa Fluor 568™	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Alexa Fluor 594™	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Atto 488	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Atto 550	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Auramine-rhodamine	SXJ-5822FLS-19002	
Azami Green	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Biosearch Blue	SXJ-5822FLS-19000	SXJ-5822FLS-39000
BODIPY FL/pH7.2	SXJ-5822FLS-19002	SXJ-5822FLS-39002
CAL Fluor® Red 590	SXJ-5822FLS-19004	SXJ-5822FLS-39004
CAL Fluor® Red 610	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Calcein	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Calcium Green™-1	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Calcofluor White	SXJ-5822FLS-19000	
Coumarin	SXJ-5822FLS-19000	SXJ-5822FLS-39000
Cy2™	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Cy3.5™	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Cy3™	SXJ-5822FLS-19004	SXJ-5822FLS-39004
DAPI	SXJ-5822FLS-19000	SXJ-5822FLS-39000
Dil	SXJ-5822FLS-19004	SXJ-5822FLS-39004
DiO	SXJ-5822FLS-19002	SXJ-5822FLS-39002
DsRed	SXJ-5822FLS-19004	SXJ-5822FLS-39004
DyLight 350	SXJ-5822FLS-19000	SXJ-5822FLS-39000
DyLight 488	SXJ-5822FLS-19002	SXJ-5822FLS-39002
DyLight 549	SXJ-5822FLS-19004	SXJ-5822FLS-39004
DyLight 594	SXJ-5822FLS-19006	SXJ-5822FLS-39010
EGFP	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Emerald GFP	SXJ-5822FLS-19002	SXJ-5822FLS-39002
FAM	SXJ-5822FLS-19002	SXJ-5822FLS-39002
FITC	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Fluo-4	SXJ-5822FLS-19002	SXJ-5822FLS-39002
FM™ 4-64	SXJ-5822FLS-19006	
FusionRed	SXJ-5822FLS-19006	SXJ-5822FLS-39010
GFP	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Hoechst 33258	SXJ-5822FLS-19000	SXJ-5822FLS-39000
JC-1 non-ratiometric	SXJ-5822FLS-19002	
Killer Red	SXJ-5822FLS-19006	SXJ-5822FLS-39010
LysoTracker Blue/MeOH	SXJ-5822FLS-19000	SXJ-5822FLS-39000
LysoTracker Green/pH 5.2	SXJ-5822FLS-19002	SXJ-5822FLS-39002
mCherry	SXJ-5822FLS-19006	SXJ-5822FLS-39010
MitoTracker Green FM/MeOH	SXJ-5822FLS-19002	SXJ-5822FLS-39002
MitoTracker Orange/MeOH	SXJ-5822FLS-19004	SXJ-5822FLS-39004
MitoTracker Red/MeOH	SXJ-5822FLS-19006	SXJ-5822FLS-39010
mKate2	SXJ-5822FLS-19006	SXJ-5822FLS-39010
mPlum	SXJ-5822FLS-19006	
mRFP1	SXJ-5822FLS-19006	SXJ-5822FLS-39010
mWasabi	SXJ-5822FLS-19002	SXJ-5822FLS-39002
Oregon Green™ 488	SXJ-5822FLS-19002	SXJ-5822FLS-39002
R-phycoerythrin	SXJ-5822FLS-19004	
ReAsH-CCPGCC	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Resorufin	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Rhod-2	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Rhodamine Red™-X	SXJ-5822FLS-19006	SXJ-5822FLS-39010
ROX	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Sulforhodamine 101	SXJ-5822FLS-19006	SXJ-5822FLS-39010
SYBR® Green I	SXJ-5822FLS-19002	SXJ-5822FLS-39002
SYTO 9/DNA	SXJ-5822FLS-19002	SXJ-5822FLS-39002
TagRFP	SXJ-5822FLS-19004	SXJ-5822FLS-39004
TAMRA	SXJ-5822FLS-19004	SXJ-5822FLS-39004
tdTomato	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Tetramethylrhodamine isothiocyanate	SXJ-5822FLS-19004	SXJ-5822FLS-39004
Texas Red®	SXJ-5822FLS-19006	SXJ-5822FLS-39010
Texas Red®-X	SXJ-5822FLS-19006	SXJ-5822FLS-39010
TRITC	SXJ-5822FLS-19004	SXJ-5822FLS-39004
TurboFP650	SXJ-5822FLS-19006	
X-rhod-1/Ca2+	SXJ-5822FLS-19006	SXJ-5822FLS-39010
ZsGreen1	SXJ-5822FLS-19002	SXJ-5822FLS-39002